

Automatic Soxhlet Fat Extractor

Model: SZF-06A/ SZF-06C

● Product Description

Fat Analyzer tests fat content according to Soxhlet Extraction principle、use weight determination method to determine the fat content。 It features extraction, reagent recovery, cooling system and temperature controlling

● Characteristic

- ▲ Small volume, fast heating speed, even heating.
- ▲ Accurate digital temperature controlling, easy to operate.
- ▲ Extraction bottle and apparatus are made of glass to avoid other spills

● Technical parameters

- ▲ Measuring range: oil content within the scope of 0.5%-60% of good, feed, oil and various kinds of fat products
- ▲ Measuring sample quantity: 6 simultaneously
- ▲ Heating time: within 10 min
- ▲ Recovery system: automatic
- ▲ Recovery rate: ≥80%
- ▲ Temperature controlling range: room temperature -100°C
- ▲ Power: AC220V/50Hz, 1000W
- ▲ Shape size: 760*230*680mm
- ▲ Weight: 25kgs



● Model difference

Model	SZF-06A	SZF-06C
Shape size	760×230×680	670×250×600
Weight	25kg	23kg
Difference	Even heating and automatic recycling	Even heating and automatic recycling, automatic soaking and extraction, save working time.