

SEED TESTING INSTRUMENT

Seed Neatness Workbench

Model: TJD-800

●Product description

TJD-800 is used come with magnifier. There are lighting devices on both the top and bottom of its observation platform, convenient for observing and cleaning seeds.

Observation platform with soft white light and even reflection can magnify testing objects, used for seeds purity testing.

●Technical parameters

Magnification: 0-7 times

Illumination power: 20-30W

Power: 220V, 50HZ

Net weight: 8KG



Seed Neatness Workbench

Model: TJD-900

●Characteristics

- ▲ Both sides of the work table adopt slope design according to the human body engineering principle; the operation is comfortable and labor-saving.
- ▲ Removable drawer is under the instrument, it is used for hold well-chosen seeds.
- ▲ Extended edge is designed on the left, right and back side of workbench to avoid seeds' lateral leakage.
- ▲ Glass surface can be drawn out to replace tube easily.
- ▲ The original seeds-fall holes are convenience for taking out the chosen seeds.
- ▲ Wooden structure, ultra thin observation platform, provided with desktop illumination magnifying glass.
- ▲ Standard configuration contains 2 different chromaticity glasses, widely use.

●Technical parameters

Magnification: 0-7 times

Surface size: 360 * 330cm

Table base size: 940 * 330cm

Luminous panel size: 300 * 270cm

