


Melting-Point Apparatus Series



WQD-400 / WQD-360 Automatic dropping-point and softening-point Apparatus



Characteristics

- HD video detection system, observing 2 samples at the same time.
- Automatically record measuring process, support video playback and data recalling.
- Compliant with ASTM D3461, ASTM D3104, ASTM D3954, ASTM D6090.
- Touch screen, user-friendly operation.
- Widely used to component control and quality test of tar , asphalt , paraffin , resin , fats , rosin , lubricant , Vaseline , medical ointment , ointment , suppository , edible grease.

Specifications

Parameters/Model	WQD-400	WQD-360
Measuring range	Room temperature-400°C	Room temperature-360°C
Measuring modes	Auto / Visual	Auto / Visual
Number of sample	2	2
Measuring methods	600	400
Number of Users	200	150
Audit trail	Optional	Optional
Resolution	0.1°C	
Linear heating rate	0.1°C / min~20°C / min	
System	Android	
Data storage	128G	64G
Display	10.1 inch touch screen, 720P HD camera	
Interface	1 x RS232 , 3 x USB , 1 x RJ-45	
Power Supply	100-240V , 50 / 60Hz , 150w	
Package Dimensions	495mm×493mm×350mm	
Weight (Gross)	5 kg	