


MHS-B-6L/MHS-B-3L Six in One Ceramic HotPlate Magnetic Stirrer



- ◆ Ceramic platform, excellent quality
- ◆ No radiation, acid-resistant, alkali-resistant, easy to clean
- ◆ Touch key, stepless speed regulation
- ◆ LED display actual temperature and time
- ◆ With time and date record, longest time: 999 minutes
- ◆ Six positions, convenient for comparing test parameters
- ◆ Widely used for heating and mixing of laboratory chemical liquids.

Model	MHS-B-3L	MHS-B-6L
Speed range	start up ~1600 rpm	start up ~1600 rpm
Heating power	3x600W	6x600W
Stirring power	3x60W	6x60W
Temperature range	50~320℃	50~320℃
Temperature precision	±0.1℃	±0.1℃
Timer	0-999min	0-999min
Working plate size(mm)	180x180mm x3	180x180mmx6
Voltage	220V 50Hz	220V 50Hz
Characteristic	Microprocessor control Digital display each temperature and time	Microprocessor control Digital display each temperature and time
Dimension(mm)	640x300x150	1288x300x150
Package(mm)	730x380x230	1370x380x230
G.W./N.W.	22/18kg	34/27kg