


TR-80 TUBER ROTATOR



- 1 Streamlined design, beautiful and elegant appearance
- 2 Long axis flip mixing method, with good uniformity and high efficiency in sample mixing
- 3 Multiple fixtures available for mixing different sample tubes
- 4 Brushless DC motor drive, high power, long service life
- 5 Low speed operation with high torque, strong power, and wide mixing range

Applications

Suitable for various experimental analyses, suspension mixing, tissue sample washing and precipitation, extraction, dissolution, etc. for comparative analysis of samples, and mixing well



Technical Parameter

Model	TR -80
Motion type	Flip over
Speed regulation range	20~80 rpm, stepless speed regulation
Overall power	20W
Work angle	/
Standard specifications	50ml*3 15ml*8 2ml*10
Power	100-240V 50-60Hz
Dimension (mm)	460×120×135
Packaging size (mm)	520x180x195
G.W/N.W	3.4/2.6kg

Long axis mixer roller:

Option 1: Centrifugal tube 50mlx12 Option 2: Centrifugal tube 15ml/10mlx23

Option 3: Centrifuge tube 1.5/2mlx28, blood collection vessel 2/5mlx28