

# CS-700



## Product Characteristics

Type	Haze Meter CS-700
Light Source	CIE-A、CIE-C、CIE-D65
Standards	ASTM D1003/D1044 ISO13468/ISO14782
Parameter	HAZE), Transmittance (T)
Spectral Response	CIE Luminosity function y
Wavelength	400-700nm
Wavelength Interval	10nm
Geometry	0/d
Measurement Area/Sample Port	16.5mm/21mm
Measurement Range	0-100%
Haze Resolution	0.01
Haze Repeatability	≤0.1
Sample Size	Thickness≤150mm
Display	5 inch TFTLCD screen
Memory	20000 values
Interface	USB
Power	220V
Working Temperature	10-40 °C
Storage Temperature	0-50°C
Size (LxWxH)	598mmX247mmX366mm
Standard Accessory	PC software (Color QC)
Optional	Measurement Fixture, Haze standard film

